



EN-E-09-001-001

## *Heating isothermal water batch*

Low price, applying in polyethylene box

Simple design, easy operation, safe application

# KTH-60T



### **MedicaTec Corporation**

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901,  
Beijing, China, 101100

TEL: + 86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicattec.cn>

E-mail: [zhangxiang@medicattec.cn](mailto:zhangxiang@medicattec.cn)



# KTH-60T

Cheap ethylene box which can be replaced regularly

Save cleansing time and trouble

Completely safe

Floating low water level detection control

Over current protection fuse

Simple operation

Button  $\Delta$  and  $\nabla$  can only be used to modify the required control temperature

Great freedom

Boxes of various sizes can be used (maximal width 460 mm)

[Applying example]

Ferment reaction

Serum thermal inactivation

Temperature maintenance for culture medium

Microorganism culture, etc.

## Product specifications and options constitute

|                                  |  |
|----------------------------------|--|
| Patten size                      | KTH-60T  |
| Water tank temperature range     | Room temperature + 5 $\square$ - + 80 $\square$                  |
| Applying environment             | 5 $\square$ - 50 $\square$ (no condense), humidity < 80%         |
| Control precision of temperature | $\pm$ 0.1 $\square$  |
| Temperature display              | Actual temperature of water tank/set temperature/digital display |
| Heating tube                     | 1000W  |
| Protective device                | Over current protection, low water level protection              |
| Tank agitation jet               | Jet efflux and agitation   |
| Occupying size                   | 316 (W) $\times$ 90 (D) $\times$ 120mm (H)                       |
| External size                    | 316 (W) $\times$ 90 (D) $\times$ 240mm (H)                       |
| Weight                           | Approx 3kg   |
| Power                            | AC220V, 7A   |

\* Spring fixing bed is optional

### MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901,  
Beijing, China, 101100

TEL: + 86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicattec.cn>

E-mail: [zhangxiang@medicatec.cn](mailto:zhangxiang@medicatec.cn)