



EN-E-09-012-001

Automatic Tube Cap Opener

◆◆◆◆◆◆◆◆◆◆ Provide the operator and the test object a safer environment ◆◆◆◆◆◆◆◆◆◆

MCO-75/MCO-100

New Release



MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901,
Beijing, China, 101100
TEL: + 86-10-89506062 Cell: +86-15110235386
URL: <http://www.medicattec.cn>
E-mail: zhangxiang@medicattec.cn

MCO-75/MCO-100 ◆◆◆◆◆◆◆◆

Easily open the tube with no manual intervention (non-glass tube)

Avoid plash contamination, and provide the operator and the test object a safer environment

【Characteristic】

- Open all tubes on the shelf simultaneously
- The caps will not contact with the opener, but be discarded into the waste bag automatically
- The opening components can be unloaded easily and washed regularly
- Opener height can be adjusted (for 6 heights)



【Specifications】

| | |
|---------------------------|---|
| TYPE | MCO-75(For 75mmTube) MCO-100(For 100mmTube) |
| Ambient | 10 ~ +50℃ (None Dewfall) , humidity less than 85% |
| Open Speed | MCO-75 15/12sec(50Hz/60Hz) MCO-100 9/7sec(50Hz/60Hz) |
| Open style | Ten tubes on the shelf can be opened simultaneously |
| Mechanical opening finger | Can be unloaded and washed (the component is exclusively designed for the tube and shelf) |
| Separator | Removable(Easy to cleaning) |
| Shooter | Removable(Easy to cleaning) |
| Dimension | 350~400 (H) ×290 (W) ×280mm (D) |
| Weight | 11Kg |
| Power | AC100V(AC220 available), 10VA |

MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901,
Beijing, China, 101100

TEL: + 86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicattec.cn>

E-mail: zhangxiang@medicattec.cn