



Microplate Jet Washer

◆◆◆◆◆◆◆◆◆◆ Used microplate Cleaning ◆◆◆◆◆◆◆◆◆◆

MW-51 II

New Release



MW-51 II



Olympus PK7300 Used microplate
Cleaning / Recycling


【Product properties】

● Cleaning / recycling speed

In the cleaning room, 3 plates for each group, total of 17 groups.
up to 51 microplates can be cleaning at one time
From washing to drying maximum 3hr for 51 plates

● 120 well microplate cleaning / recycling

Washing liquid, water and DI water
Warm wind drying



【Specifications】

TYPE	MW-51 II
Dimension	1120 (W) × 750 (D) × 1050mm (H)
Weight	250Kg
Power	AC100V, 2500VA

MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901,
Beijing, China, 101100

TEL: + 86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicattec.cn>

E-mail: zhangxiang@medicatec.cn