



EN-E-10-003-001

ELISA 96ch Microplate washer with automatic dredge function



Three set of nozzle



MW-96CR

New Release



MW-96CR

Automatic dredge function to clean nozzle

【Characteristic】

- 96ch washing at same time
- Automatic dredge function to clean the nozzle
- Automatic Head Rinsing
- Three set of nozzle, reduce contamination transfer between well

- Automatic Piston Plunging System:
After washing action, automatically cleaning washing head, to prevent fluid caused crystallization have in washing nozzle, [Piston Plunging Unit] take the initiative to clean up the foreign body in nozzle. [Piston Plunging Unit] working at the certain time on washing process, to clean up the foreign body in nozzle.



96ch 3 sets nozzle

【Specifications】

TYPE	MW-96CR
Ambient	10 ~ +50℃ (None Dewfall) , humidity less than 85%
Dispense volume	50μL~990μL
Wash time	1-9time
Wash method	96ch dispense at same time, suction at same time
Microplate	All sorts of 96ch microplate
Residual liquid	Every well≤1.5μL
Washing liquid dunk time	maximum 999 second
Washing mode	Normal, Dispense, Aspirate, Dispense End, 9 user program memory
Dripping precision	CV<6% (@100μL Dispense)
Dimension	440 (W) ×490 (D) ×260mm (H)
Weight	25Kg
Power	AC100V(AC220 available), 100VA

MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901, Beijing, China, 101100

TEL: + 86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicathec.cn>

E-mail: zhangxiang@medicatec.cn